

Downstream Processing of Antibodies

Introduction

A biological product, such as a protein or antibody, is purified and recovered by a series of processes known as "downstream processing" after it has been created by a living thing or a cell culture. The goal of downstream processing is to remove and purify the target product from other substances like contaminants, host cell proteins, and DNA that are present in the crude mixture. Several unit activities, such as clarifying, purification, polishing, and formulation, are frequently involved in downstream processing. The medicinal efficacy and market value of the final product obtained through downstream processing can be considerably impacted by its quality and yield. The successful manufacture of high-quality biological products, including various kinds of antibodies, depends on optimizing downstream processing.

Because of their intricate structure and large molecular weight, antibodies present particular difficulties during downstream processing. Impurities such host cell proteins, DNA, and aggregates are a significant problem since they can have an impact on the end product's effectiveness, safety, and stability. This can be avoided by using a combination of different purification methods to eliminate these contaminants, including affinity chromatography, ion exchange chromatography, and size exclusion chromatography. Another difficulty is the potential loss of antibody activity during downstream processing, which can be reduced by utilizing gentle purification procedures and improving the composition of the buffer. The heterogeneity of antibody glycosylation can also be problematic, although this can be overcome by creating more effective and targeted glycan removal techniques.

Downstream Methods & How They Are Applied

The separation, purification, and concentration of target molecules like proteins, enzymes, and antibiotics from a complex mixture of biologically produced components are steps in the bioprocessing process known as downstream processing. For downstream procedures, several techniques are frequently used, including filtration, chromatography, centrifugation, and precipitation.

Using a porous membrane or filter, filtration entails separating the target molecule from contaminants based on variations in surface charge or size exclusion. For this, ultrafiltration and microfiltration are frequently employed. Although microfiltration separates molecules based on particle size, ultrafiltration separates them based on their molecular weight. Both methods involve applying pressure to a solution as it passes through a filter, maintaining the target molecule while filtering out contaminants.

Another popular technique for downstream processing is chromatography. It is a method that divides mixture components according to their chemical and physical characteristics. This technique can separate target molecules with high purity and is very selective. Depending on the characteristics of the target molecule, many chromatography methods can be applied, including ion exchange, affinity, size exclusion, and hydrophobic interaction.

A method for separating molecules according to their density and size is centrifugation. When a combination is spun at a high speed, the lighter components stay at the top while the denser components fall to the bottom of the tube or centrifuge. With this method, it is possible to separate the target molecule from big particles and cellular waste.

Finally, precipitation entails the addition of a reagent to the solution to cause the target molecule to precipitate out of solution. Proteins and nucleic acids are frequently initially purified using this technique.

Advantages & Advancements

The purification of antibodies was a difficult and drawn-out procedure before the development and improvements of downstream processing. Precipitation, dialysis, and centrifugation were early procedures that could only produce a little amount of yield and purification. This frequently led to impure antibody products that had low potency and ineffective effectiveness. The manufacture of large quantities of antibodies was further constrained by the absence of effective downstream processing techniques, making it challenging to satiate the rising demand for these molecules in both research and clinical contexts. Higher purity and yield, shorter manufacturing times, and more scalabilities have been made possible by the development of downstream processing techniques including chromatography and membrane filtering. These advancements have significantly enhanced the development and commercialization of therapeutic antibodies, leading to the production of a wide range of highly effective and life-saving antibody-based drugs. The main advantages include: Higher purity, improved potency, higher yield, increased scalability, reduced production time, and enhanced comparability.

Disadvantages & Concerns

Although downstream processing for antibodies has several benefits, there are still obstacles and restrictions that must be overcome. The possibility of losing antibody function during the purification process, which can be brought on by elements like harsh pH or temperature conditions, extended exposure to air, or the presence of denaturants, is a significant problem. For

therapeutic purposes, this may affect the final antibody product's efficacy and potency. The heterogeneity in the glycosylation patterns of antibodies, which might impact their pharmacokinetic and pharmacodynamic characteristics, is another difficulty. Moreover, the manufacture and marketing of some antibodies may be constrained by the expense and complexity of downstream processing. Furthermore, the yield of the target antibody can be low due to the inefficiencies of the purification process, which can increase the cost of production. The downstream processing also requires specialized equipment and trained personnel, which can add to the cost and complexity of the process.

The need for more research and development to improve antibody downstream processing methods is generally highlighted by these problems. Increases in purification techniques, such as the use of gentler conditions or new purification techniques, can lessen activity loss and boost yield. To guarantee constant product quality, improvements in analytical techniques can also enable better monitoring and control of the downstream process. Moreover, the adoption of alternative expression platforms, such as those based on plant or insect cells, may open up fresh opportunities for increasing the effectiveness and cost-effectiveness of antibody downstream processing.

Possible Disadvantage Mitigations

To reduce the problems related to antibody downstream processing, many options and modifications are being researched. One strategy is to monitor the purifying process in real-time using more sophisticated analytical techniques, which would allow for improved process control and the early detection of any faults. Another tactic is to purify the antibody under milder conditions that are less likely to denature or harm the molecule. Other expression methods are

also being researched, including cell-free systems, yeast, and algae, which may have benefits for yield, scalability, and affordability. Moreover, improvements in glycoengineering technology might make it possible to produce antibodies with more dependable glycosylation patterns, which might enhance the antibody's pharmacokinetic and pharmacodynamic qualities. The use of machine learning and artificial intelligence to optimize downstream processes is a current topic of research that has the potential to increase yield, lower costs, and reduce unpredictability.

To reduce the problems and enhance the downstream processing for antibodies, several modifications and options are being investigated. The creation of new purification methods that are more effective, less harmful, and capable of catching a wider variety of antibody types is one potential option. This can entail the application of brand-new resins, membranes, or other substances as well as brand-new elution and buffer exchange techniques. Moreover, efforts are being made to enhance the effectiveness and reproducibility of downstream processing through developments in automation, robotics, and process control technologies. The utilization of continuous processing techniques is also gaining popularity since it has the potential to increase yields, shorten processing times, and reduce product variability. Finally, efforts are being made to create more thorough analytical techniques for characterizing and monitoring antibody products throughout the downstream processing pipeline, which may aid in early problem detection and purification process optimization for each individual antibody product.

By identifying the ideal composition and concentration of nutrients needed for the growth and productivity of the cell line, medium optimization is utilized to enhance this process. Using both statistical experimental designs and empirical techniques, this is accomplished. The statistical designs allow for the screening of various parameters and their interactions, which may be used to forecast the best formulation. Examples of these statistical designs include response

surface approach and design of experiments. To ascertain how different nutrient concentrations, affect cell growth and productivity, the empirical approaches evaluate how well cells operate in media with various nutrient concentrations. High productivity, quick cell development, and reliable product quality should all be supported by the ideal medium. Once the ideal medium has been found, it can be further adjusted using process management techniques including media supplementation and feeding schemes to boost output and product quality.

Conclusion

In conclusion, over the past ten years, the productivity of monoclonal antibody manufacturing has significantly increased due to improvements in upstream processes, making mAb production more practical and affordable for therapeutic applications. This development has been facilitated by the optimization of host cells, bioreactor technology, and process parameters, the creation of new bioreactor designs and the ongoing investigation into.

Some challenges that are faced regarding upstream processing of mAbs, is the optimization of cell culture conditions to ensure high productivity and product quality, while maintaining cell viability and preventing contamination. This requires a lot of work in understanding the cellular physiology of cell production. Also, looking at the use of advanced techniques such as metabolic flux analysis and omics approaches [5]. Another challenge is the scalability of the process. Upstream processing typically involves the use of large-scale bioreactors, which can be difficult to control and maintain at high cell densities [4]. This requires careful process optimization and monitoring, as well as the use of advanced control strategies such as fed batch or perfusion culture.

References

- [1] Jain, E., & Kumar, A. (2008). Upstream processes in antibody production: Evaluation of Critical Parameters. *Biotechnology Advances*, 26(1), 46–72.
<https://doi.org/10.1016/j.biotechadv.2007.09.004>
- [2] Duvar, S., Hecht, V., Finger, J., Gullans, M., & Ziehr, H. (2013). Developing an upstream process for a monoclonal antibody including medium optimization. *BMC Proceedings*, 7(S6).
<https://doi.org/10.1186/1753-6561-7-s6-p34>
- [3] Wahlgreen, M. R., Meyer, K., Ritschel, T. K. S., Engsig-Karup, A. P., Gernaey, K. V., & Jørgensen, J. B. (2022). Modeling and simulation of upstream and downstream processes for monoclonal antibody production. *IFAC-PapersOnLine*, 55(7), 685–690.
<https://doi.org/10.1016/j.ifacol.2022.07.523>
- [4] Yang, J.-D., Lu, C., Stasny, B., Henley, J., Guinto, W., Gonzalez, C., Gleason, J., Fung, M., Collopy, B., Benjamino, M., Gangi, J., Hanson, M., & Ille, E. (2007). Fed-batch bioreactor process scale-up from 3-L to 2,500-L scale for monoclonal antibody production from Cell Culture. *Biotechnology and Bioengineering*, 98(1), 141–154. <https://doi.org/10.1002/bit.21413>
- [5] Kurz, M. (n.d.). Upstream processing: Regulatory considerations regarding quality aspects of monoclonal antibodies. *BioPharm International*. Retrieved April 7, 2023, from <https://www.biopharminternational.com/view/upstream-processing-regulatory-considerations-regarding-quality-aspects-monoclonal-antibodies>